


ANCHORMAN™

CORTICAL TIBIAL FIXATION



Product Information

ANCHORMAN - CORTICAL TIBIAL FIXATION

ANCHORMAN is a tibial fixation device for Anterior Cruciate Ligament (ACL) reconstructive surgery.

BENEFITS OF ANCHORMAN

STRONG CORTICAL FIXATION

ANCHORMAN provides exceptionally strong cortical fixation. Its fixed loop PEEK design offers up to 98% less cyclic displacement compared to both traditional tibial cortical fixation devices and interference screws^{1,2}.


SIMPLE TO USE. ONE SIZE

ANCHORMAN is simple to use. One size fits all tibial tunnels between 6-10mm in diameter and only requires, one simple step for tibial tunnel preparation.

PRESERVATION OF THE ACL TIBIAL REMNANT

ANCHORMAN allows preservation of the native tibial fibres. With a fixed PEEK loop, ANCHORMAN offers a less traumatic fixation with no threads or moving sutures to potentially damage the graft.

CYCLIC DISPLACEMENT (MM)



- ANCHORMAN
- ARTHREX TIGHTROPE RT
- S&N BIOSURE PK



PRODUCT CODES

ANCHORMAN™	
Code	Description
D201019	ANCHORMAN Tibial Hook Implant
D201039	Tibial Template

Ref 1 Aga, C, Biomechanical Comparison of Interference Screws and Combination Screw and Sheath Devices for soft Tissue Anterior Cruciate Ligament Reconstruction on the Tibial Side, American Journal of Sports Medicine, 2013
 Ref 2 Robert, H, A comparison of four tibial fixation systems in hamstring, graft anterior ligament reconstruction, European Journal of Orthopaedic Surgery & Traumatology, 2014

